Biotest HYCON

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Biotest AG – Landsteinerstrasse 5-
63303 Dreieich - Germany

Biotest HYCON

RCS Compressed Gas Adapter
(CG-Adapter) for RCS Airborne Sampler

Art. – No.: 940 720

Contents		Page
0.	Scope of delivery	10
1.	Application	11
2.	Principles of Operation and Construction	12
3.	Operating Instructions	12
	3.1. General Notes	12
	3.2. Sterilization and Decontamination before Use	12
	4.3. Preparation of RCS Gas Sampler	13
	4.4. Performance of Compressed Air Sampling	13
4.	Accessories	13 - 14
5.	Technical Data	14
6.	Warranty	14

Scope of delivery

- 1x Compressed gas adapter
- 1x Manual
- 1x Transport case

1. Application


The Biotest RCS Compressed Gas Adapter in conjunction with the RCS Plus, RCS High Flow, and RCS Isolator Microbial Air Samplers, is used to investigate the microbiological content of pressurized gases according to ISO14698. It is routinely used in areas where specified air microbial counts must not be exceeded due to legal requirements or individual standards.

The RCS Compressed Gas Adapter (RCS CG-Adapter) allows for a simple yet effective measurement of both high or low concentrations of air microbes in the following environments:

- **Isolators and laminar flow**
- **Sterile and aseptic production lines in pharmaceutical industry**

2. Principles of Operation and Construction

The RCS Compressed Gas Adapter consists of a diffusor, reduction nozzle and a tube connection, where the

compressed gas pipeline can be connected to the adapter.

The tube connection on the RCS CG-Adapter has a thread size of R1/4 and is designed to provide a universal connection.

If the tube connection needs to be replaced by a different size or a quick coupler connection, make sure that it is securely attached to the reduction nozzle - an important component of the RCS CG-Adapter.

Please note that during any maintenance on the device all pieces must be maintained in a safe environment and should not be lost.

As the gas flow enters the RCS CG-Adapter, through the tube connection, the velocity of the volume flow is restricted by the reduction nozzle. The conduit of the air stream then expands in the diffusor allowing the gas to arrive inside the attached RCS Microbial Air Sampler almost pressure free. Excess air from the compressed gas line flows over the protection cap of the RCS Microbial Air Sampler through the small gap between the diffusor and protection cap.

Information on the operation and sampling principles of the RCS Microbial Air Sampler can be found in the Operations Manual of the respective unit.

According to ISO14698 a sample volume of 1 m³ must be used for microbial counts in gasses.

3. Operating instructions

3.1. General Notes

Standard equipment for evaluating the microbial content of compressed gasses includes the RCS CG-Adapter, an RCS Microbial Air Sampler (RCS Plus, RCS High Flow and RCS Isolator), one or more calibrated rotors and agar strips.

Due to the dry characteristic of compressed gasses, it is recommended that TCI-γ agar strips be used for this evaluation.

In order to avoid contamination of your sample, sterile gloves should be worn when handling the RCS CG-Adapter and the microbial air sampler. Please ensure that all air outlet openings of both units are not obstructed during operation in order to avoid turbulences or variations in flow rate.

It is also recommended that the agar strips be placed inside the respective rotors and in the RCS airborne sampler in a sterile environment such as a laminar flow hood.

3.2. Sterilization and Decontamination before Use

An occasional thorough cleaning of the RCS CG-Adapter is recommended. Sterilization and decontamination measures should be established for the adapter and routinely performed. Since the adapter is made of aluminum it can be rinsed or sprayed with a detergent or commercially available disinfectant with a neutral pH (pH= 6.5 – 7.5). The neutral pH prevents harm to the anodic coating of the aluminum. Alternatively, the unit can be autoclaved with the rotor and protection cap of the RCS Microbial Air Sampler at 121°C. After sterilization it is recommended that the Compressed Gas Adapter be stored in a sterile container or bag.

3.3. Preparation of the RCS Air Sampler

Prior to microbial sampling, the RCS Air Sampler must be assembled with the compressed air adapter and an equipped rotor in a sterile environment. This can be accomplished by gently pushing the adapter over the protection cap of the RCS and allowing it to be secured in this position by the clamps located on the RCS CG-Adapter. Should a number of tests be performed in one sampling location, it is recommended to equip a corresponding number of rotors with agar strips and transport them in a sterile container to the sampling location. According to ISO 14698 a sample volume of 1m³ must be used for evaluating the microbial content of gases. Program your microbial air sampler to sample 1000 litres and select the proper calibration factor of the used rotor.

3.4 Performance of Compressed Air Sampling

Before sampling the pressurized gas line make sure that you have assembled the RCS CG-Adapter and the RCS Microbial Air Sampler in a sterile environment and then connect the adapter to the compressed air pipeline. The instrument must then be placed horizontally on a solid surface on its tripod. Please ensure a tight and close fit between the adapter and compressed air pipeline. If necessary, disinfect the connecting parts. Next open the pressure valve making sure that it does not exceed 1 bar + 10%. As the air begins to flow through the adapter a marked hissing noise will come from the expansion cylinder.

The gas should escape for approximately 1 minute, before switch on the RCS Air Sampler (see paragraph 3.3.) After completion of air sampling, close the compressed air valve and remove RCS Air Sampler from expansion cylinder.

Transfer instrument or rotor with agarstrip in a sterile container to a sterile environment. Remove agar strip there,

slide back into strip cover, label and process for incubation. If multiple microbial counts are to be performed at the same compressed air outlet, the compressed air valve should not be closed during rotor replacement.

4. Accessories

For microbial counts in compressed gases, a RCS Air Sampler is needed.

4.1. RCS Air Sampler

4.1.1. RCS Plus

(Order No. 940 310): High performance air sampler for sampling volumes up to 1m³ with a flow rate of 50 ltr./min. incl. accessories in a carrying case.

4.1.2. RCS High Flow

(Order No. 940 215): High performance air sampler for sampling volumes up to 1m³ with a flow rate of 100 ltr./min. incl. accessories in carrying case

4.1.3. RCS Isolator

(Order No. 940 810): High performance air sampler with detached operating and sampling unit for sampling volumes up to 1m³ and a flow rate of 100 ltr./min. incl. accessories in carrying case.

4.2. Agar Strips

4.2.1. Agar Strips TCI-γ

(Order No. 941 125); γ-irradiated, double wrapped casein soybean digest agar with neutralizers for determination of total counts in aseptic environments and also in peroxid-containing air. Also suitable for microbial counts with dry gases.

5. Technical Data

Diffusor and Reduction Nozzle: made from anodized aluminum. Autoclavable and disinfectable with commercially available detergents, H₂O₂-sterilizable connection parts: stainless steel.

Range of application: 1 bar \pm 10%. Higher pressures must be reduced by pressure regulators. Reduction nozzles for lower pressures available upon request. Suitable for non-flammable and non-toxic gases.

Flow rate: With a primary pressure of 1 bar, the flow rate is approx. 180 ltr./min.

Sample volume: the recommended amount of air for air sampling is 1 m³ according to ISO14698.

6. Warranty

Biotest declares its instruments to be free from defects of materials and workmanship. The period of warranty is one (1) year commencing on the date of commissioning of the instrument as evidenced by means of an invoice or other document.

The warranty covers the rectification of all faults which occur as a proven direct result of faulty materials or workmanship during the warranty period.

The following items are not covered by the warranty:

1. Parts subject to normal wear and tear.
2. Malfunction caused by accidental or intentional damage, or the use of parts from other manufacturers, or through repairs or modifications performed by persons not authorized by Biotest to carry out such work.
3. Defects which occur as a result of normal wear and tear or from using the equipment under conditions not specified in the instructions for use.

The user is not entitled to any remedy in law towards Biotest, howsoever arising, apart from those remedies mentioned here. Biotest accepts no liability or warranty towards the user. The legal rights of the purchaser towards the seller remain unaffected by this.

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Biotest Italia S.r.l.
Biotest (Schweiz) AG
Biotest Diagnostic Corp.