Biotest MMB Prep Salmonella

Manual

Version 06.09

Biotest MMB Prep Salmonella (2x 50 Preps)
Art. No. 944 505

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hycon@biotest.de

Introduction

Salmonellosis is a food-borne infection caused by several strains of the bacterial species *Salmonella enterica* (serovar). The consumption of contaminated food derived from animals or plants is the main reason for Salmonellosis.

The real-time PCR method, in contrast to classic microbiological culture methods, enables the rapid and sensitive testing of raw material and processed food for absence of *Salmonella*.

The **Biotest MMB Prep Salmonella** kit can be used for template preparation for PCR analysis.

Intended Use

The kit is intended to be used for the isolation of Salmonella DNA from enrichment cultures, suspensions, swabs and food samples. The isolated DNA is suitable for use with the Biotest MMB Salmonella real-time PCR test.

Kit Components (50 preparations) and Storage

- 1x Lysis Buffer (20 ml)
- 1x Binding Buffer (10 ml)
- 1x Wash Buffer (60 ml, after adding 42 ml ethanol (absolute, not included in the kit))
- 1x Elution Buffer (10 ml)

- 1x Receiver Tube 2 ml (50 x)
- 1x Receiver Tube 1.5 ml (50 x)
- 1x Spin Filter (50 x)

All kit components are stored dry and at room temperature (14-25 °C).

Additional required reagents

Ethanol (absolute, >96 %)

Additional required equipment and disposables

- Microcentrifuge
- Thermomixer / heating block (100 °C)
- Adjustable pipettes
- Sterile DNA-free filter tips
- 1.5 ml and 2 ml microcentrifuge tubes

Quality Control

According to the Quality Management System, each lot of Biotest MMB Prep Kits is tested against predetermined specifications to ensure consistent product quality.

Warning and Precautions

Not for *in vitro* diagnostic use.

- Caution:** The test must be performed by well-trained and authorized laboratory technicians.
- Caution:** All reagents should be handled in accordance with good laboratory practice taking appropriate precautions.
- Caution:** Biological samples can contain infectious agents. Follow all appropriate precautions, e.g., the regular check of gloves for contamination with sample material.
- Caution:** Do not use reagents past the expiration date printed on the label.
- Caution:** Do not use reagents with any evidence of turbidity or microbial contamination.

For more information, please consult the appropriate material safety data sheets (MSDS). These are available online in PDF format at www.biotest.de. In the section Microbiological Tests & Systems you can find, view, and print the MSDS for each Biotest MMB Prep kit and kit component.

Protocol

Before using the kit for the first time add 42 ml ethanol absolute to the Wash Buffer and mix thoroughly.

Before each isolation procedure the required amount of Elution Buffer is transferred to a reaction tube and pre-warmed to 50-60 °C.

Preparation of starting material

1. Transfer 1 ml of the enrichment culture or corresponding specimen under sterile conditions into a 1.5 ml reaction tube (not provided with the kit).
2. Centrifuge for 5 min at 14,000 x g.
3. Remove and discard supernatant by pipetting.

Lysis of starting material

4. Add 400 µl* of **Lysis Buffer**, mix briefly and close the reaction tube tightly.
5. Incubate under continuous shaking for 10 min at 99 °C.

Setting of optimal binding conditions

6. Centrifuge for 1 min at 14,000 x g.
7. Transfer approx. 300 µl of the liquid supernatant into a new 1.5 ml reaction tube (not supplied with the kit).

* The volume depends on the specific sample matrix. It is possible to vary the volume up to 1.0 ml. For further questions please send contact to hycon@biotest.de.

Binding of the DNA to a Spin Filter

8. Add 200 µl **Binding Buffer** to the filtrate and mix the sample.
9. Place a **Spin Filter** into a yellow **2 ml Receiver Tube**.
10. Transfer the solution directly onto the filter and incubate at room temperature for 1 min.
11. Centrifuge for 1 min at 14,000 x g.
12. Discard the flow-through and place the spin filter back into the same receiver tube.

Purification of the bound DNA

13. Add 550 µl **Wash Buffer** to the spin filter and centrifuge for 1 min at 14,000 x g. Discard the flow-through and place the spin filter back into the same receiver tube.
14. Add 550 µl **Wash Buffer** to the spin filter and centrifuge for 1 min at 14,000 x g. Discard the flow-through and place the spin filter back into the receiver tube.

Drying of the Spin Filter

15. Remove the residual ethanol by final centrifugation for 2 min at 14,000 x g.

Elution of nucleic acids from the Spin Filter

16. Place the spin filter into a clear **1.5 ml Receiver Tube**.
17. Add 100 µl of the preheated **Elution Buffer** (50-60 °C) directly onto the spin filter.
18. Incubate 3 min at room temperature.
19. Centrifuge for 1 min at 10,000 x g. After centrifugation discard the spin filter.

The eluted DNA is ready-to-use for PCR analysis. The DNA can be stored for up to 24 hours at 4 °C. For a storage time of more than 24 hours it should be kept at -20 °C.

Additional Information

For any inquiries or to place an order, contact your local Biotest distributor.

Domestic customers in Germany may contact Customer Service at:

Phone: +49 (0) 6103 801-496

Fax: +49 (0) 6103 801-505

mail@biotest.de

For technical inquiries please contact our Technical Support at:

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Biotest (Schweiz) AG
Biotest Italia S.r.l.
Biotest Diagnostics Corporation